




## VELA Samba 510



Pos	Pcs	Beskrivelse	Description	Beschreibung	Stock no
1	1	Gasfjeder chrome yderrør 300 N 200 mm	Gascolumn 300N, 200 mm, chrome	Gasfedersäulen 300 N, 200 mm, chrom	091030200
2	1	1 bæreplade, HS, åben rotation	1 seat support, Rs, open rotation	1 trägerplatte, Ls,	18012130999999
3	1	Beslag for rygbøjle S510	Bracket for backrest stem S510	Klammer für rückenstäbe S510	926180
4	1	Ryghængsel T50	Backrest hinge T50	Rück scharnier T50	925617
5	1	Ryg for S510 Bore R.7	Back for S510 Bore R.7	Rück für S510 Bore R.7	08147511672
6	1	Rygbøjle for Samba 510	Backrest stem S510	Rückenstäbe S510	07381
7	2	Afstandskive Ø20x8	Distance disc Ø20x8	Distanzscheibe Ø20x8	924886
8	1	Håndskrue VCT.40 p-M8X16	Handscrew VCT.40 p-M8X16	Handschraube VCT.40 p-M8X16	925553
9	2	Skrue med invendig sekskant CH sort krom	Screw M6 x 16 black	Schraube M6 x 16 Schwarz	8001106016
10	1	Sæde for taburet Ø34 cm Bore R1.2 f/ fast ryg	Seat for stool Ø34 cm Bore R1.2 f/ fast ryg	Sitz für hocker Ø34 cm Bore R1.2 f/ fast ryg	06432420672
11	6	DIN 912 - M6 x 20	DIN 912 - M6 x 20v	DIN 912 - M6 x 20	8001206020
12	1	Rørben 15x30 mm flad-oval	Plug 15x35 mm	Zwinge 15x35 mm	922125
13	1	Knopfjeder for 10 mm rør	Pop-out pin for 10 mm tube	Druck-stifte Ø10 mm Rohr	923362
14	6	Skive	Washer	Scheibe	8802207
15	1	Stel m/hjul Alu fod typr 1, Blank alu, N/	Frame with wheels	Untergestell mit Rollen	020129906
15.1	1	Alu fod type 1N blank	Frame with out wheels	Untergestell ohne Rollen	030129999
15.2	5	Hjul Dobb. 2, bløde m/ 11x20, 50 mm sort	Wheel soft 50 mm black	Rollen weich 50 mm Schwarz	012202